Invia® Motion™

NEGATIVE PRESSURE WOUND THERAPY SYSTEM

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Assistance

Please read the entire Instructions for Use before trying to operate this device.

If you have any questions about your Invia Motion Negative Pressure Wound Therapy (NPWT) System, please contact your healthcare professional at the number below:

Healthcare professional contact information:

Please keep the Patient Instruction for Use in an easily accessible location.

Table of Contents

EN Instructions for Use	3
ES Instrucciones para uso	42
FR Mode d'emploi	81
General Warnings.....	4
Purpose of the Invia NPWT Pump Device.....	7
Intended Use.....	7
When Not to Use Device (Contraindications).....	8
Service Life.....	8
Description of the Device.....	9
Display.....	10
Conditions that Affect Use.....	11
General Safety at Home.....	12
Sleeping.....	13
Showering, Bathing and Washing.....	14
When You Leave the House.....	14
Set Up Instructions.....	15
Set Up Carrying Case.....	17
Charge Battery.....	19
Change Canister/Tubing Set.....	21
Notifications and Troubleshooting.....	24
Info Used Pump Time.....	31
Cleaning/Disinfection.....	33
Warranty.....	33
Maintenance.....	33
Disposal.....	33
Technical Specifications.....	34
Signs and Symbols.....	35
Technical Data.....	37
Electromagnetic Compatibility.....	37

General Warnings

WARNING

The individual responsible for operating this device on a daily basis (patient, caregiver) must be able to read, understand and follow the directions provided in this Instructions for Use manual, on the pump display and by the healthcare professional. The individual must be able to hear and respond to the auditory signals.

Only use the Invia Motion NPWT System on the person for whom it was ordered and only for its intended use.

Do not perform the therapy without your healthcare provider's supervision.

WARNING

Visually check the wound dressing and the Invia Motion NPWT Pump frequently.

If you see large amounts of blood in the dressing, tubing or canister:

1. Turn off the pump by pressing the  button and holding for 3 seconds.
2. Apply pressure to your wound.
3. Seek immediate emergency medical assistance.

If you have difficulty breathing, seek immediate emergency medical assistance.

 WARNING

Immediately contact your healthcare professional if you experience any of the following symptoms:

- Significant vision or hearing loss
- Pain, soreness or redness around the wound
- Swelling, itching or rash around the wound
- Fever
- Strong odor from the wound
- Thick green drainage or any change in color of drainage
- Nausea, vomiting or diarrhea
- Difficulty breathing
- Dizziness or confusion
- Allergic reactions like redness, rash or hives, swelling or severe itching

 WARNING

Once the healthcare professional has set the pump pressure and mode settings, do not press the pressure  or mode  setting buttons.

Pressing these buttons may change the settings of the pump and could cause harm.

WARNING

The pump must be used 24 hours per day. Do not stop the pump unless instructed by your healthcare professional or if excessive bleeding occurs.

If the pump is stopped for more than the time frames shown below, the dressing must be changed. Contact your healthcare professional to request a dressing change.

Gauze dressings: Change the dressing if the pump is stopped for more than 24 hours.

Foam dressings: Change the dressing if the pump is stopped for more than 2 hours.

WARNING

If the pump stops working and cannot be switched on, immediately call your healthcare professional.

CAUTION: Federal law restricts this device to sale by or on the order of a physician.

Purpose of the Device

Your healthcare provider has prescribed Negative Pressure Wound Therapy (NPWT) for you.

Negative Pressure Wound Therapy has been shown to help promote healing in several different kinds of wounds.

Negative pressure (suction) is delivered to a wound through a pump. The suction, or negative pressure, helps to remove fluid and close your wound. Your healthcare professional will place a special dressing on your wound and a tube will go from your wound to the canister on the pump. The pump turned on and set to the prescribed pressure to remove the fluid from your wound. The fluid will be collected in the canister. The treatment will be applied 24 hours a day.

Intended Use

The portable Medela Invia Motion Negative Pressure Wound Therapy (NPWT) System is indicated to create an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudates and infectious material. It is intended for use in hospitals, clinics, home care (HC) and long term care (LTC) settings on adult patients with chronic, acute, subacute, traumatic, dehiscent wounds, partial thickness burns, ulcers (such as diabetic, neuropathic, pressure or venous insufficiency), flaps and grafts.

When Not to Use Device (Contraindications)

The Invia Motion NPWT System is contraindicated in the presence of:

- Malignancy in the wound
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Necrotic tissue with eschar present
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

If you have concerns, speak with your healthcare professional.

WARNING

Therapy to be performed by a healthcare professional only.

Service Life

The service life of the Invia Motion Endure is three years.

Description of the Device

Check the delivery package of the Invia Motion NPWT system for completeness and general condition.

WARNING

Do not use any other accessories than the ones described below.

- a 1 Invia Motion NPWT Pump
- b 1 Invia Motion Canister/Tubing Set 0.15l (150 ml) (sterile, single use)
- c 1 Invia Motion Power Supply
- d 1 Invia Motion Carrying Case
- e 1 Invia Motion Negative Pressure Wound Therapy System Clinician Instructions for Use
- f 1 Invia Motion Negative Pressure Wound Therapy System Patient Instructions for Use

Display

Buttons to Operate the Invia Motion NPWT Pump

Do Not Press
(only for healthcare professionals)

Do Not Press
(only for healthcare professionals)

The image shows the front panel of the Invia Motion NPWT Pump. It features a central LCD display. Above the display are two buttons: 'Change Pressure Setting' (with a bell icon) and 'On/Off' (with a power icon). Below the display are two buttons: 'Constant/Intermittent Mode' (with a waveform icon) and 'Mute Acoustic Signal/Info Used Pump Time' (with a speaker icon crossed out). The Medela logo is at the bottom.

On/Off

Mute Acoustic Signal/
Info Used Pump Time (Not applicable for REF 087.4004 Invia Motion Endure)

The image shows the back panel of the Invia Motion NPWT Pump. It features a large label with technical specifications, a barcode, and a 'Run-time: yy days' indicator. An arrow points from the 'Display Information' table to the top of the pump.

Display Information

Leakage	Canister Full
Battery Status	Info Used Pump Time
Pump Is Running/ On-State In Intermittent Mode	Off-State In Intermittent Mode

Conditions that Affect Use

CAUTION

The Invia Motion NPWT Pump must remain in an upright position or lying display side up during use. Do not invert the pump or turn the display side down. Please see below.

General Safety at Home

CAUTIONS / WARNINGS

- Invia Motion NPWT Pump does not detect blockage along tubing or connections. Check Invia Motion tubing frequently throughout the day to make sure tubing is not kinked or blocked. If blockage is identified along Invia Motion tubing, switch Invia Motion NPWT Pump off and contact your healthcare professional.
- Do not use extension cords with the Invia Motion NPWT Pump.
- Position the power cord so that it will not cause someone to trip and fall.
- Never place the power supply cord around your neck.
- Keep the power supply cord away from hot surfaces.
- Do not unplug the power supply cord from the wall outlet by pulling on the cord.
- Do not try to fix or repair the Invia Motion NPWT Pump. Contact your caregiver if the pump is not working.
- Keep the Invia Motion NPWT Pump clean and dry.
- Never place the Invia Motion NPWT Pump in water or liquids.
- If the Invia Motion NPWT Pump gets wet, dry with a towel. Do not dry the Invia Motion NPWT Pump in a microwave.
- Keep the Invia Motion NPWT Pump away from children and pets.
- Keep the Invia Motion NPWT Pump protected from direct sunlight.

- The use of mobile telephones, local area networks including wireless, walkie-talkies (two-way radios), and cordless telephones can affect the Invia Motion NPWT Pump. This may cause insufficient battery charging or may lead to an internal fault. See the chapter Notifications and Troubleshooting for actions to take in the case of an internal fault.

Certain mobile phones may interfere with the Invia Motion NPWT system when receiving and making calls close to the pump. Prior to using a mobile phone near the Invia Motion NPWT system, check that the pump continues to function while making or receiving a call.

If you observe a change in the vacuum level or an error condition, move the Invia Motion NPWT system at least 3.3 feet (1 meter) from the mobile phone when making or receiving calls. Turn the Invia Motion NPWT system on and off again to reset it after an error condition.

Whenever possible, a safety distance of at least 3.3 feet (1 meter) between the Invia Motion NPWT Pump and a mobile phone is recommended.

Sleeping

- Make sure the Invia Motion tubing will not become kinked.
- Prevent the Invia Motion NPWT Pump from falling to the floor during sleep.
- Always charge the Invia Motion NPWT Pump while sleeping.

Showering, Bathing and Washing

- Do not use the Invia Motion NPWT Pump while bathing or showering.
- The Invia Motion NPWT Pump can be disconnected for a short period of time (up to 20 minutes) in order to change clothes or perform personal hygiene. Discuss bathing plans with your healthcare professional.
- Prevent the Invia Motion NPWT Pump from getting wet.
- Do not place the Invia Motion NPWT Pump where it can fall into a shower, sink or tub.
- Do not touch the Invia Motion NPWT Pump if it has fallen into water. Unplug the device immediately from the electrical outlet and disconnect the pump from the dressing.

When You Leave the House

- Check the battery status and recharge the battery if necessary before leaving the house.
- For your convenience, place the Invia Motion NPWT Pump into the carrying case and store excess tubing safely in the tubing pouch located in the case.
- Do not forget to take a spare Canister/Tubing Set and the power supply with you. The power supply can be stored in the tubing pouch.

Set Up Instructions

WARNING

Check the Invia Motion NPWT Pump frequently to make sure it is running. Visually check the wound dressing, the area around the wound and the fluid in the canister for signs of infection or other complications as provided by your healthcare professional.

Switch Invia Motion NPWT Pump on

Press

Display view when the pump is running. Arrow is blinking ►.

WARNING

Do not try to perform therapy without the supervision of your healthcare professional.

WARNING

NPWT must run 24 hours a day. Do not stop the therapy unless instructed by your healthcare professional or if excessive bleeding occurs.

WARNING

The dressing should be changed and therapy restarted by the healthcare professional if therapy is discontinued for:

- More than two 2 hours when using a foam dressing.
- More than 24 hours when using a gauze dressing.

Switch Invia Motion NPWT Pump off

Press the  button and hold for 3 seconds. Press only once.

Backlight

The backlight can be activated by pressing the  button. Without further interaction the backlight will turn off after 10 seconds.

Set Up Carrying Case

The Invia Motion NPWT Pump is intended to be worn in a carrying case during operation. One is included in the product package.

1. Open the carrying case.
2. Slide the Invia Motion NPWT Pump into the designated upper pocket and close the Velcro strap.
3. Place the excess tubing into the designated lower pocket and close the Velcro strap.
4. Fold the pump and tubing sections of the carrying case together and align the inner Velcro tabs to secure.
5. Secure the bag by closing the lid and fastening the button.

The carrying case strap can be worn over the shoulder or around the waist.

To wear over the shoulder, attach the strap clips to the carrying case. To wear as a belt, thread the strap through the two loops located on the back of the case.

Shoulder strap

Belt strap

Charge Battery

When the battery is fully charged, it should last approximately 10 hours. If the Invia Motion NPWT Pump is running continuously due to an air leakage, the battery may need to be recharged earlier (after approximately 4 hours).

Battery Status Symbols

full

half full

fully discharged (empty)

If the battery becomes fully discharged (empty), please follow the instructions below.

Open the Direct Current (DC) socket by pulling the small protective cover open. Connect the small end of the power supply into the Invia Motion NPWT Pump. Plug the large end of the power supply into a wall outlet. The moving bars on the battery display will indicate that charging is occurring.

The battery will be fully charged after approximately 4 hours. Once charged, the  symbol will appear constantly on the display. After charging is completed, unplug the power supply from the wall and disconnect the Invia Motion NPWT Pump from the power supply.

Close the Direct Current (DC) socket by pressing the small protective cover over the opening.


CAUTION

NPWT therapy must be continued while charging the pump's battery. Do not turn the pump off while charging the battery.

Always charge the battery overnight.

Direct Current (DC)
socket

Change Canister/Tubing Set


 To be used by healthcare professionals or adequately trained caregivers only.

Change the canister when it appears full or when the canister full signal sounds or when the canister full symbol  is displayed.

1. Wash hands.
2. Switch the Invia Motion NPWT Pump off by pressing the  button and holding for 3 seconds. Press only once.

3. Disconnect the Invia Motion tubing connector from the wound tubing connector by separating the pump from the patient tubing.

Have a paper towel or tissue available to cover the tip of the used tubing to protect against wound fluid leakage.

4. Remove the pump from the carrying case and remove the full Canister/Tubing Set by rotating the canister to the left.

5. Wash hands.
6. Obtain a new Invia Motion Canister/Tubing Set and open the packaging.


CAUTION

Do not use an Invia Motion Canister/Tubing Set if the sterile packaging is damaged.

7. Remove the new Canister/Tubing Set from the packaging and attach it to the Invia Motion NPWT Pump by pressing the canister to the pump and rotating to the right.

8. Connect the Invia Motion Tubing Connector to the wound drain tubing.

9. Switch the Invia Motion NPWT Pump on by pressing the  button. Slide the pump into the carrying case.

10. Discard the used Canister/Tubing Set according to instructions provided by your healthcare professional.

Notifications and Troubleshooting

Do not panic if an acoustic signal sounds.
Follow the instructions below.



1. If the Invia Motion NPWT Pump is in the carrying case open the lid of your carrying case.

2. Press the  button to mute the acoustic signal for 60 seconds (it can be pressed every time the signal sounds, with the exception of "Battery Empty" and after the countdown of 15 or 60 days plus 72 hours of safety time is used up).




3 Check the display for one of the following visual symbols.

Visual Symbol displayed	What does this symbol mean?	What do I have to do?
 Blinking! + acoustic signal!	The canister is full	Mute the acoustic signal by pressing the  button.   1. Wash hands. 2. Switch the Invia Motion NPWT Pump off by pressing the  button for 3 seconds. Press only once. 3. Disconnect the Invia Motion Tubing Connector from the wound tubing connector by separating the pump from the patient tubing. Have a paper towel or tissue available to cover the tip of the used tubing to protect against wound fluid leakage. 4. Remove the pump from the carrying case. Remove the full Canister/Tubing Set by rotating the canister to the left. Discard the used Canister/Tubing Set according to instructions provided by your healthcare professional. 5. Wash hands.

Visual Symbol displayed	What does this symbol mean?	What do I have to do?
 Blinking!  acoustic signal!	The canister is full	6. Obtain and open a new Invia Motion Canister/Tubing Set package.  Caution Do not use an Invia Motion Canister/Tubing Set if the sterile packaging is damaged. 7. Remove the Canister/Tubing Set from the packaging and attach it to Invia Motion NPWT Pump. 8. Connect the Invia Motion Tubing Connector to the wound drain tubing. 9. Switch the Invia Motion NPWT Pump on by pressing the  button. Slide the pump into the carrying case. 10. Ask your clinician how to properly dispose of a full canister. If you are unable to change the canister contact your healthcare professional.


Visual Symbol displayed	What does this symbol mean?	What do I have to do?
 Blinking! + acoustic signal!	There is a leak present	Mute the acoustic signal by pressing the  button. Check the dressing and tubing connections for leaks and pat the transparent film (clear bandage) placed over the wound and surrounding skin.  If you are unable to eliminate the leak, switch the Invia Motion NPWT Pump off by pressing the  button and holding for 3 seconds. Contact your health-care professional for instructions.

Visual Symbol displayed	What does this symbol mean?	What do I have to do?
 Blinking! + acoustic signal!	The battery is low	 Mute the acoustic signal by pressing the  button. The battery requires charging. Connect the Invia Motion NPWT Pump to the supplied electrical plug and power source. Moving bars on the battery display will indicate that charging is occurring. Do not turn the pump off while charging.  If you do not have access to electricity and the battery is fully discharged (empty), switch the pump off by pressing the  button.

Visual Symbol displayed

What does this symbol mean?


What do I have to do?


+
acoustic
signal!

The battery is fully discharged (empty)



The pump will only work for ten more minutes.

Charge the battery immediately. Connect the Invia Motion NPWT Pump to the supplied electrical plug and power source. Moving bars on the battery display will indicate that charging is occurring. Do not turn the pump off while charging.

If you do not have access to electricity and the battery is fully discharged (empty), contact your healthcare professional for instructions.

Visual Symbol displayed	What does this symbol mean?	What do I have to do?
60 Blinking! (Example of 60 Day version)	Available operation time elapsed	Contact your healthcare professional. Do not turn the pump off! 
E + acoustic signal!	Internal fault	Turn the Invia Motion NPWT Pump off by pressing the  button and holding for 3 seconds.  Contact your healthcare professional immediately for further instructions.

If you are unable to identify any symbol on the display, switch off the pump by pressing the  button and holding for 3 seconds. Contact your healthcare professional.

WARNING

If the pump fails or cannot be switched on, call your healthcare professional immediately.

Info Used Pump Time

NE

REF 087.4006 Invia Motion – 15 Day Negative Pressure Wound Therapy System


REF 087.4000 Invia Motion – 60 Day Negative Pressure Wound Therapy System

These two Invia Motion NPWT Pump versions have a run-time of 15 and 60 days.

Press the  button to obtain information regarding the treatment time used.

The picture below illustrates a 60 days run-time Invia Motion pump that has operated for 6 days and 23 hours. The remaining run-time of this Invia Motion NPWT Pump would be 53 days and 1 hour.

Please inform your healthcare professional when the remaining operation time is three days.

REF 087.4004 Invia Motion – Endure Negative Pressure Wound Therapy System

This Invia Motion NPWT Pump version has no countdown.

Cleaning/Disinfection

NE

The Invia Motion NPWT Pump can be wiped with a damp cloth using a mild household cleaner.

Warranty

The warranty period is for 2 years after the date of delivery, but for 15 day and 60 day versions expires after 15 or 60 days of usage plus 72 hours of safety time, if used in accordance with Instructions for Use. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage or use by unauthorized persons.

Maintenance

The Invia Motion NPWT Pump is maintenance free and will not require service.

If an Invia Motion NPWT Pump fails within the warranty period due to a manufacturing defect, it will be replaced. The original pump will need to be returned to the supplier.

Disposal

Ask your healthcare professional how to dispose of used Canister/Tubing Sets and the Invia Motion NPWT Pump at the end of your therapy.

Technical Specification

Accessories Storage Conditions

Accessories must remain in the packaging and be stored within a temperature range of -4° Fahrenheit to +122° Fahrenheit (-20° Celsius to +50° Celsius).

Operation Environmental Conditions

The Invia Motion NPWT Pump and accessories must be operated within a temperature range of +41° Fahrenheit to +104° Fahrenheit (+5° Celsius to +40° Celsius). Do not operate the product in extreme cold or heat. Avoid high relative humidity (15–93 % is permissible).

CAUTION

Do not operate the Invia Motion NPWT Pump above 9,842 feet above sea level (3000 meter above sea level, unless in the pressurized passenger cabin of an aircraft).

Protection Class IP22

The Invia Motion NPWT Pump is protected against ingress of solids larger than .50 inches (12.5 millimeter) and dripping water.

Signs and Symbols

This Symbol indicates to consult the Instructions for Use.

This Symbol indicates a CAUTION or WARNING associated with the Device (see page 3).

This Symbol indicates a class II device.

This Symbol indicates not to dispose of the Device with Household Refuse (for EU only).

This Symbol indicates the Date of Manufacture (four Digits for the Year and two Digits for the Month).

This Symbol indicates the Name and the Address of the Manufacturer.

This Symbol indicates the Device is sterilized using ethylene oxide.

This Symbol indicates the device is MR Unsafe.

This Symbol indicates that this is a Prescription Device. CAUTION: Federal Law restricts this Device to sale by or on the order of a licensed healthcare practitioner. (for US only).

This Symbol indicates a type BF applied part.

This Symbol indicates the manufacturer's Catalog Number.

This Symbol indicates the manufacturer's Serial Number.

This Symbol indicates the manufacturer's Batch Code.

This Symbol indicates the device has Protection against the ingress of dripping water and solid foreign objects of .50 inches and greater.

This Symbol indicates that the Device should not be used after the End of the Year and Month shown.

This Symbol indicates to follow Instructions for Use.

This Symbol indicates the C TUV US NRTL marking of the Device (equivalent to UL and CSA Mark).

This Symbol indicates not to use the Device if package is damaged.

This Symbol indicates the number of items that the package contains.

pcs

This Symbol indicates the number of items.

This Symbol indicates the Direct Current Socket.

This symbol indicates a Single Use Device. Do not reuse the Device.

This Symbol indicates the Temperature Limitation for Operation, Transport and Storage.

This Symbol indicates the Humidity Limitation for Operation, Transport and Storage.

This Symbol indicates the atmospheric Pressure Limitation for Operation, Transport and Storage.

This Symbol indicates to keep the Device dry.

This Symbol indicates to handle the fragile Device with care.

This Symbol indicates to keep the Device away from sunlight.

This Symbol indicates a carton package.

This Symbol indicates that the material is part of a recovery/recycling process.

⁰¹²³ This Symbol indicates that the Device is in Conformance with the Medical Device Directive 93/42/EEC (for EU only).

Technical Data

Contact your local representative or Medela if you have questions regarding the technical information below.

Electromagnetic Compatibility

IEC 60601-1-2, Table 1

EMC

Electromagnetic Emissions		
The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT System should assure that they are used in such environment. Contact your local representative or Medela if you have questions regarding the technical information below.		
Emission Tests	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Invia Motion NPWT System uses RF energy only for their internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment. The Invia Motion NPWT System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

WARNING

The Invia Motion NPWT System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Motion NPWT System should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Immunity

The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. The customer or the user of Invia Motion NPWT System should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 16.5 kV air without power supply ± 8 kV air with power supply	The relative humidity should be at least 5 %.
NOTE: At high levels of ESD, the system may reset itself. The Invia Motion NPWT Pump will restart on its own at correct settings without interruption of therapy and requiring no action from the user.			
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for input / output lines not applicable	Mains power quality should be that of a typical domestic or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode not applicable	Mains power quality should be that of a typical domestic or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	$<5\% U_T$ ($>95\%$ dip in U_T) for 0.5 cycle $40\% U_T$ (60% dip in U_T) for 5 cycles $70\% U_T$ (30% dip in U_T) for 25 cycles $<5\% U_T$ ($>95\%$ dip in U_T) for 5 sec	Mains power quality should be that of a typical domestic or hospital environment. If the user of the Invia Motion NPWT System requires continued operation during power mains interruptions, it is recommended the Invia Motion NPWT System be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Separation distance $r = \frac{I}{188}$ Where I is the current in amperes in a power bus or an appliance wire and r is the separation distance in meters.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Electromagnetic Immunity

The Invia Motion NPWT System is intended for use in the electromagnetic environment specified below. Users of the Invia Motion NPWT System should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Invia Motion NPWT System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0.4 \sqrt{P}$
Radiated RF IEC 61000-3	3 V/m 80 MHz to 2.5 GHz	With power supply: 3 V/m 80 MHz to 800 MHz 3 V/m 800 MHz to 2.5 GHz Without power supply (battery): 10 V/m 80 MHz to 800 MHz 10V/m 800 MHz to 2.5 GHz	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Invia Motion NPWT System are used exceeds the applicable RF compliance level above, Invia Motion NPWT System should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Invia Motion NPWT System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and The Invia Motion NPWT System The Invia Motion NPWT System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Users of the Invia Motion NPWT System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Invia Motion NPWT System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	m		
W	150 kHz to 80 MHz $d = 0.4 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.04	0.12	0.23
0.1	0.13	0.38	0.73
1	0.4	1.2	2.3
10	1.3	3.8	7.3
100	4	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Medela and Invia are registered trademarks of Medela Holding AG.
Avance is a registered trademark of Mölnlycke.
Medela AG
Medela Inc. USA
Medela Canada Inc.