

INVACARE®
ROLLITE®
ROLLATORS



The **Invacare® Rollite® Rollator** line is unlike any other rollator line in the industry. Invacare has taken the Rollite Rollator one step further, offering the sizes, wheel types and colors that consumers want. The Rollite line offers adult and junior sizes in different colors. For improved mobility over rougher terrain, large 8" wheel models are available in adult and tall sizes (Midnight Blue color). Along with additional features, the Rollite line provides the comfort, durability and quality that consumers have come to expect from Invacare products.

65100R Rollator Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

HCPCS Code: E0143/E0156




1. Folds into a compact unit for storage. Built-in seat handle allows for easy transport



65100 Rollator Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

HCPCS Code: E0143/E0156



1. Folds into a compact unit for storage. Built-in seat handle allows for easy transport



NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

68100 Rollator Features

- Lightweight aluminum frame
- Flip-up seat with built-in handle
- Ergonomic dual-paddle folding mechanism
- Folds into a compact unit for storage and transport
- Flexible backrest for comfort
- Ergonomic handbrakes lock the rear wheels for security
- Adjustable seat-to-floor height

HCPCS Code: E0143/E0156



1. Folds into a compact unit for storage. Built-in seat handle allows for easy transport



1. **Basket**

- Model # 65110 for 65100, 65100R and 68100-TA
- Model # 65109 for 65100-JR, 65100R-JR and 68100
- 10 lb. weight capacity
- Position basket underneath seat to secure valuables



2. **Tray**

- Model # 65408
- For food and beverage placement
- For use on all Rollite models

NOTE: Rollators are not intended to be used as wheelchairs.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further details.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
©. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-234
111258



Model #	Dimensions Overall	Dimensions Inside Handgrip	Dimensions At Base	Depth Folded	Seat Dimensions	Seat-to-Floor Height	Wheels
65100, 65100R	27" W	18.25" W	22"Wx24"D	7"	17"Wx9.5"D	19"-23"	5"
65100-JR,65100R-JR	27" W	18.25" W	22"Wx24"D	7"	17"Wx9.5"D	16.5"-20.5"	5"
68100	27" W	18.25" W	22"Wx26"D	7"	17"Wx9.5"D	19.5"-23.5"	8"
68100-TA	27" W	18.25" W	22"Wx26"D	7"	17"Wx9.5"D	22"-26"	8"



	Height Adjustments	Weight	Weight Capacity	Client Height*	Quantity Per Carton	Limited Warranty
65100, 65100R	32"-36"	13.5 lb.	300 lb.	5'4"-6'	1	5 Years Frame 1 Year Components**
65100-JR,65100R-JR	29"-33"	13.5 lb.	300 lb.	4'11"-5'7"	1	5 Years Frame 1 Year Components**
68100	32"-36"	13.5 lb.	300 lb.	5'4"-6'	1	5 Years Frame 1 Year Components**
68100-TA	35"-39"	16.5 lb.	300 lb.	5'11"-6'6"	1	5 Years Frame 1 Year Components**

* Client height is approximation. Each individual should always consult with his/her physician or therapist to determine proper adjustment.

** No tires/wheels

NOTE: Rollators are not intended to be used as wheelchairs