

Spenco® Gel Insoles

Cat. No. 39-818



Target Consumer:

Consumers seeking all day comfort for their feet. Consumers standing for long periods of time: workers, students and travelers.

Recommended For:

Daily wear – absorbing shock to reduce foot fatigue and sore feet.

For workers, acts like an anti-fatigue mat in the shoe.



Best For Use In:

Any casual or athletic shoe.



Components*:

TPR Gel (triple density),
Heel thickness – 8.35mm [3/8"],
Forefoot thickness – 4.31mm [5/32"],



Instructions:

Remove existing insole (some brands of shoes glue the insole to the shoe, but careful removal will not damage the shoe). Insert the Spenco® Gel Insole into the shoe. Trim with scissors if necessary.



Care:

Hand wash with warm water. Air dry.

Features & Benefits:

Triple Density Design comforts, cushions and supports – advanced design for ultimate comfort.

- 1 **TPR Comfort Gel** thermoplastic rubber combines superior cushioning and energy return for ultimate comfort.
- 2 **Metatarsal Arch Support** relieves pressure from the ball-of-foot.
- 3 **Top Cloth with Ultra-Fresh Silpure Antimicrobial** helps eliminate odor causing bacteria and reduces friction.
- 4 **Stability Cradle** supports the arch and heel for greater comfort. Makes it easier to place into the shoe.
- 5 **Extra Heel Cushioning** softer gel in the heel absorbs impact at heel strike.

One-Year Unconditional Guarantee

Sizing Chart:

Size	0	1	2	3	4	5	6
U.S. Women's	3/4	5/6	7/8	9/10	11/12		
U.S. Men's			6/7	8/9	10/11	12/13	14/15
U.K.	2/3	4/5	5/6	7/8	9/10	11/12	13/14
European	32/34	36/38	38/40	40/42	42/44	44/46	46/48

Call toll-free: 1-800-877-3626
In Canada: 1-800-387-9538

 www.spenco.com



*Measurements are approximate and vary by size.